# FERMILAB/NICADD PHOTOINJECTOR LABORATORY



H. Julia Brito SIST Program Tuesday, August 6, 2002

#### PERSONAL BACKGROUND



- Undergraduate at The University of Texas in Austin
- Mechanical Engineering
- Concentration in Biomedical Engineering


#### TEMP HOME @ FERMI

- Location
  - A- Zero
  - Beams Division
  - Research and Development lab
- Family
  - DESY = Deutsches Elektronen-Synchrotron (Germany)
  - NICADD = Northern Illinois Center for Accelerator Detector and Development
  - LBNL = Lawrence Berkeley National Laboratory
  - University of Rochester

#### **OUTLINE**

- Brief tour of A0 Photoinjector beam line
- Short discussion of work with diagnostic apparatus
- Concise summary of future goals
  - Photoinjector lab
  - Personal

## PHOTOINJECTOR BEAM LINE SCHEMATIC




#### KLYSTRON FOR PHOTO-CATHODE GUN




- Output power = 3 MW
- Located outside of beam line cave
- Radio Frequency (RF)
   waves @ 1.3 GHz travel
   via a waveguide—
   minimal losses for
   higher frequencies

#### ENTERING THE CAVE



- •RF waves produce electric field in 1.5 cell cavity
- Phase timing is crucial for electrons produced by cathode prep chamber to be accelerated
- Solenoids encircle the cathode gun for electron compression and focus
- Something needed to initially pick off electrons

#### **CLASS IV LASER**




- Laser beam is created in oscillator cavity
- Many components, more things to go wrong
- Problems are not-soobvious (lamp strike failure)



#### INFRARED TO ULTRAVIOLET LIGHT

Multi-pass stages used to amplify energy of the laser







IR is converted to UV via doubling crystals and is guided into the cave

#### LASER ENTERS THE BEAM LINE


Quantum Efficiency = # of electrons / photon



Laser optics in cave (Remote-control Iris)



Bad phase scan



Good phase scan


#### NINE CELL SUPER-CONDUCTING CAVITY



\*Each cell progressively steps up the acceleration of the electrons


\*Energy of the electrons increases from ~4 MeV to ~17 MeV

#### CRYOGENIC MONITOR OF 9 - CELL




#### **MAGNETS**

- <u>Dipole magnets</u>: used for steering
- <u>Quadrupole magnets</u>: used for focusing; focuses in the transverse plane i.e. focuses in dimension of space
- <u>Chicane</u>: consists of 4 dipole magnets; steers beam down, straight, up, & straight again; focuses in the longitudinal plane i.e. focuses in dimension of time



#### END OF THE BEAM LINE

Spectrometer = magnet used to measure the energy at the end of the beam line



Total Energy = 
$$sqrt((1.54*I_{spec})^2 + 0.511^2)$$

_MULT_	3.000	MULT CONTROLLER	į	3.000		N/A	
SOLMN	1.000	SOLENOID MAIN	*	179.999	s	A	
SOLBUC	1.000	SOLENOID BUCKING	*	189.998	s	A	
SOLSEC	1.000	SOLENOID SECONDAR!	*	180.003	s	A	
PSPECT	0.200	SPECTROMETER PS	*	9.496	s	A	
CAMCTL	1.000	AVMS CNTL	*	lundefined	s		EN
FCUPX2	1.000	ACT#1 FCUPX2	, <b>*</b>	lundefined	s		OUT
FLAGX2	1.000	ACT#2 FLAGX2	*	lundefined	s		out
FLAGX3	1.000	ACT#3 FLAGX3	*	lundefined	s		OUT

# **TURNING ON!**



#### **DIAGNOSTICS**

- A-Zero Control Room
  - Parameter pages
  - Oscilloscope readings
  - Data acquisition
- Cameras measure transverse bunch profile
- Interferometer measures longitudinal bunch length



#### INTERFEROMETER @ XL4

Best approach would be to work without windows i.e. have the interferometer be a part of the beam line

However this type of interferometer would:

\*be hard to align


\*cost ~\$50,000

#### Solution:

Insert flag in beam line @ 45° angle and design a customized viewport

## VIEWP

- Crystalline quartz window
- Stainless steel recessed tube and flange
- Vacuum tight conditions






#### INTERFEROMETER DESIGN

Designer:
Uwe Happek

@ The Univ.
of Georgia



Purpose:
detectors use an
inverse Fourier
transform to reflect
shape of electron
distribution in bunch


#### TRANSLATION STAGE




- Controller/keypad/motor stage system
- Communicate with manufacturers for correct setup
- Labview interface

STABLE 1


- Drawing
- Machining
- Polishing
- Packaging
- Shipping





#### **FUTURE GOALS...**

#### ...of A-Zero Photoinjector Lab



#### ...of a little Texan intern



Biomedical Optics, Cell/Tissue Engineering, Nanotechnology, Biomedical Instrumentation, Prosthetic Engineering, etc.