The background features abstract, colorful swirls in shades of green, purple, and blue, interspersed with small yellow triangles. The text is centered and has a subtle drop shadow.

# Multi-Straw Prototype Detector *Tension Measurements*

*Oswaldo A. Lozoya*  
*University of Texas at El Paso*



# Straw Chamber

▶ The Straw Chamber is a detector whose wires pass through single straws, connected to a high voltage source. Inside the straw, two incoming wires are joined by a non-conductive glass fiber bead.



# Straw Chamber cont'd

- Straw (cathode): Kapton - Al - Kapton (LR) layers; 4 mm diameter
- Wire (anode): Au-coated, 25- $\mu\text{m}$  W wire



# Straw Chamber vs. Wire Chamber

- Straw chamber has been chosen for the BTeV project for one main reason: functionality
  - Broken wires





# Straw Prototypes

- Two kinds: Single and Multi-Straw
- Certain tension in the wires is required for the detector to operate properly
- Prototypes have a single continuous wire, e.g. no glass fiber bead
- Resonant frequency  $\propto$  Tension

# How to make a prototype



# How to make a prototype cont'd



# Theory

- In a vibrational approach, the wire represents a continuous system (many DOF's)
- Due to its characteristics (tensioned and fixed), max amplitude  $\propto$  max power occurs at ***harmonic frequencies***
  - ***Wave interference and superposition***
  - ***Nodes and antinodes***



# Theory cont'd

Governing equation:

$$f = \frac{n}{2L} \sqrt{\frac{T}{\mu}}$$

$$n = 1, 2, 3, \dots$$

***f*** = harmonic frequency (Hz)

***T*** = tension (N)  
***2<sup>nd</sup>...***

***n*** = harmonic (1<sup>st</sup>,

***L*** = length of wire (m)

***μ*** = linear density of wire (kg/m)



# Tension tests


- Function Generator (**FG**)

- An AC current is induced in the wire, located in a magnetic field
- When the wire is driven to resonance, the effect is perceived in an oscilloscope
- Current flow induction necessary (**not possible in real detector** due to wire discontinuity and bead insulation)

- Tone Generator (**TG**)

- The wire is connected to a constant potential source
- The straw-wire system resembles a cylindrical coaxial capacitor;  $Q=CV$
- The wire is vibrated with a loudspeaker, and the power output is read and displayed in a computer (LabVIEW VI) for different frequencies
- No current needs to be induced; capacitance changes suffice

# Tension tests cont'd



*Function Generator (FG)*



*Tone Generator (TG)*

# TG Hardware / Software


- Hardware

- MSP / SSP
- Voltage source ( $\sim 81$  V)
- Amplifier
- Loud Speaker
- AC signal circuit
- Oscilloscope
  - Tektronix TDS 350
  - Two channel, 200 MHz, 1 Gs/s


- Software

- Tone Gen 1.1 (LabVIEW)
  - Generates tone
  - Plots Power as function of frequency
    - Resonance = Max Power

# TG Hardware




# TG Software: VI Panel






# TG Software: VI Diagram




# TG Setup





# TG Setup cont'd






# Experiments

- TG vs. FG
  - Evaluate the difference between FG and TG readings
  - 34 samples
  - MSP used for such purpose
  - ***Different tensions mean different frequencies***
- TG accuracy
  - Evaluate how consistent and accurate is TG
  - 20 samples
  - SSP used in the study
  - ***Expected frequency known (approx. 250 Hz)***


# Results: TG vs. FG



- Mean = 1.80 Hz
- Std. Dev. = 1.23
- Graph tends slightly more to a small difference, as shown
- Good news: both of them can be trusted in prototypes *if they are right*
- ***Fact: only TG can be used in real detectors!!***

# Results: TG accuracy

- Mean = 250.66 Hz
- Expected  $\approx$  250 Hz
- Std. Dev. = 0.42
- ***Good news: TG is working right!!***
  - ***TG is potentially useful for real detector***





# Conclusion

- TG is a trustworthy device for wire tension measurement, that can be used in the detector and prototypes
- The process has been speeded up, but must be faster if possible for real detector (80,000+ wires)